



ENERTECH

Electro Surgical Generator

Ensurg SSE - TUR Plus

Enertech providing it's widest range of high performance state of the art COST EFFECTIVE Electro Surgical Generators available today in the market. The machines come with Instant Tissue response technology and Tissue impedance measurement systems. These models offer top of the line features like 99 program modes, multiple blend modes, endo cut facility, best monopolar cutting and coagulation performance with precision controlled bipolar output and host of safety features with adherence to EC 60601 guidelines. The range also includes low cost analog models for price conscious customers.



Ensurg SSE - TUR Plus

Ensurg SSE - TUR is a 400W analog model with number of cut, blend and coagulation modes with to do under water surgery. This is the model with all around capabilities and simple of operation. It is the lowest cost model in similar class available today in the market.

Performance Features :-

➤ Patient Plate Return Electrode Monitoring System :

Patient Plate Return Electrode Monitoring System provides instant indication of the quality of contact of the patient plate return electrode. Failure or poor contact results in fault indication with sound alarm which substantially reduces the risk of burn under the patient return electrode.

➤ High Power Efficiency of 97%

The power efficiency of Enertech generators are the highest in it's class 97%. This leads to multiple benefits like increased reliability, most compact size, high portability, economy and reduced power consumption.

Technical Specification

	Max.Power(Watts)	RatedLoad(Ohms)
Operating Frequency	480 KHz	-
Pure Cut	400 W	500
Blend 1	250 W	500
Blend 2	200 W	500
Spray	120 W	500
Fulguration	150 W	500
Desiccation	150 W	500
Bipolar Coagulation	80 W	100
Dimension (W x H x D)	300 mm x 150 mm x 340 mm	
Weight	5.5 Kg.	
Input Voltage	170 to 260 VAC 50 Hz	

Standard Accessories

Monopolar Double Pedal Footswitch
/ Bipolar Single Pedal Footswitch



Patient Plate With Cable Cord
Autoclavable



Silicon Patient Plate (Optional)



Chuk Handle With Electrode
Autoclavable



Bipolar Forcep Autoclavable
With Cable



Disposable Hand Switch Pencil
Imported



Ensurg SSE - TUR Plus

