



Endoscopy Cameras

EMCAM



EMCAM

Performance Features :-

- Completely medical degree design, soaked in a liquid to sterilize.
- Camera head has 2.5m Over Length cable for easy operation.
- Support hand-actuated, automatic white balance.
- Excellent optics adaptor ,breathable, water leakage, remove mist, rotate.
- Adjusting window size to adapt to different thickness of the Endosope.
- Cheap, cost-effective.
- Cater to the sophisticated needs.
- Brightness adjustment.

Technical Specification

| Model | EMCAM |
|-------------------|---------------------|
| Camera Head | 1/3" 1 Chip. |
| Resolution | PAL |
| Scanning Standard | 480 Lines. |
| Outputs | BNC |
| SNR | 60dB |
| Shutter | Automatic |
| White Balance | Automatic / Manual. |

EMCAM

